



PS-Based Fluorescent Liposome, Dextran-pHrodo Green Encapsulated

Catalog: Lipo-042RG

PRODUCT INFORMATION

Name	PS-Based Fluorescent Liposome, Dextran-pHrodo Green Encapsulated
Cat No.	Lipo-042RG
Product Overview	<p>pHrodo Green is a pH sensitive dye which increases in fluorescence as pH decreases from more neutral to acidic. Dextran-pHrodo Green conjugated liposomes are commonly used for detecting endocytosis and phagocytosis, respectively, as the dyes display little to no fluorescence outside the cell at neutral pH and upon internalization, emit a bright red or green signal due to the acidic environment of endosomes and phagosomes. Excitation/Emission (nm): 509/533</p>
Lipid Composition	<p>DOPS/DOPC/CHOL(2.0/5.0/3.0 mM) DOPS: 1,2-dioleoyl-sn-glycero-3-phospho-L-serine (sodium salt)/CAS: 90693-88-2 DOPC: 1,2-dioleoyl-sn-glycero-3-phosphocholine/CAS: 4235-95-4 CHOL: Cholesterol/CAS: 57-88-5</p>
Form	Liquid
Storage Buffer	20 mM HEPES, pH 7.4
Concentration	Lipid Concentration 10mM
Stability	2 months
Storage	2 °C to 8 °C (do not freeze)