



# POPE/CHOL/POPG Liposomes, Carboxyfluorescein Encapsulated

**Catalog: Lipo-058RG**

## PRODUCT INFORMATION

---

<b>Name</b>	POPE/CHOL/POPG Liposomes, Carboxyfluorescein Encapsulated
<b>Cat No.</b>	Lipo-058RG
<b>Product Overview</b>	Carboxyfluorescein is a self-quench dye and liposomes which encapsulate carboxyfluorescein can be used for fluorescent leakage assay. This liposome is made to order and delivered with no dialyzed and cleaned up. Customers should prepare the proper spin column for cleaning up the prep. It is advised to use the products within 30 days of the manufacturing date. Excitation/Emission (nm): 492/517 nm
<b>Lipid Composition</b>	POPE/CHOL/POPG (3.34/1.34/0.33 mM) POPE: 1-palmitoyl-2-oleoyl-sn-glycero-3-phosphoethanolamine/CAS:268550-95-4 CHOL: Cholesterol/CAS: 57-88-5 POPG: 1-palmitoyl-2-oleoyl-sn-glycero-3-phospho-(1'-rac-glycerol) (sodium salt)/CAS: 268550-95-4
<b>Form</b>	Liquid
<b>Storage Buffer</b>	20 mM HEPES, pH 7.4
<b>Concentration</b>	Lipid Concentration 5mM
<b>Stability</b>	1 months
<b>Storage</b>	2 °C to 8 °C (do not freeze)

---