



DOPC/CHOL Liposomes, Carboxyfluorescein Encapsulated

Catalog: Lipo-064RG

PRODUCT INFORMATION

Name	DOPC/CHOL Liposomes, Carboxyfluorescein Encapsulated
Cat No.	Lipo-064RG
Product Overview	Carboxyfluorescein is a self-quench dye and liposomes which encapsulate carboxyfluorescein can be used for fluorescent leakage assay. This liposome is made to order and delivered with no dialyzed and cleaned up. Customers should prepare the proper spin column for cleaning up the prep. It is advised to use the products within 30 days of the manufacturing date. Excitation/Emission (nm): 492/518 nm
Lipid Composition	DOPC/CHOL(70/30 molar ratio) DOPC: 1,2-dioleoyl-sn-glycero-3-phosphocholine/CAS: 4235-95-4 CHOL: Cholesterol/CAS: 57-88-5
Form	Liquid
Storage Buffer	20 mM HEPES, pH 7.4
Concentration	Lipid Concentration 5mM
Stability	1 months
Storage	2 °C to 8 °C (do not freeze)
