



# Negatively PG-Based Fluorescent Liposome, Dil Labeled

**Catalog: Lipo-110RG**

## PRODUCT INFORMATION

---

<b>Name</b>	Negatively PG-Based Fluorescent Liposome, Dil Labeled
<b>Cat No.</b>	Lipo-110RG
<b>Product Overview</b>	Dil Labeled Fluorescent Liposome incorporate a lipophilic dye (Dil) inside their membranes. This product could be used for tracking and detection. Excitation/Emission (nm): 549/565
<b>Lipid Composition</b>	18:1 (?9-Cis) PG/DOPC/CHOL(2.0/5.0/3.0 mM) 18:1 (?9-Cis) PG: 1,2-dioleoyl-sn-glycero-3-phospho-(1'-rac-glycerol) (sodium salt)/CAS: 67254-28-8 DOPC: 1,2-dioleoyl-sn-glycero-3-phosphocholine/CAS: 4235-95-4 CHOL: Cholesterol/CAS: 57-88-5
<b>Form</b>	Liquid
<b>Storage Buffer</b>	PBS, pH 7.4
<b>Concentration</b>	Lipid Concentration 10mM
<b>Stability</b>	2 months
<b>Storage</b>	2 °C to 8 °C (do not freeze)

---