



# Negatively PS-Based Fluorescent Liposome, Dansyl lipid (headgroup) Labeled

**Catalog: Lipo-126RG**

## PRODUCT INFORMATION

---

<b>Name</b>	Negatively PS-Based Fluorescent Liposome, Dansyl lipid (headgroup) Labeled
<b>Cat No.</b>	Lipo-126RG
<b>Product Overview</b>	18:1 Dansyl PE has the Fluorescent group at the head group of the lipid. This product could be used for tracking and detection. Excitation/Emission (nm): 336/513
<b>Lipid Composition</b>	L-alpha-PS/L-alpha-PC/CHOL(2/3/5 mM) L-alpha-PS: L-alpha-phosphatidylserine/CAS: 383908-63-2 L-alpha-PC: L-alpha-phosphatidylcholine/CAS: 97281-47-5 CHOL: Cholesterol/CAS: 57-88-5
<b>Form</b>	Liquid
<b>Storage Buffer</b>	PBS, pH 7.4
<b>Concentration</b>	Lipid Concentration 10mM
<b>Stability</b>	2 months
<b>Storage</b>	2 °C to 8 °C (do not freeze)

---