



## Positively DOTAP-Based Fluorescent Liposome, Dansy l lipid (headgroup) Labeled

Catalog: Lipo-152RG

## PRODUCT INFORMATION

Name	Positively DOTAP-Based Fluorescent Liposome, Dansyl lipid (headgroup) Labeled
Cat No.	Lipo-152RG
Product Overview	18:1 Dansyl PE has the Fluorescent group at the head group of the lipid. This product could be used for trackin g and detection. Excitation/Emission (nm): 336/513
Lipid Composition	DOTAP/DOPC/CHOL(0.5/6.5/3 mM) DOTAP: 1,2-dioleoyl-3-trimethylammonium-propane (chloride salt)/C AS: 132172-61-3 DOPC: 1,2-dioleoyl-sn-glycero-3-phosphocholine/CAS: 4235-95-4 CHOL: Cholesterol/CA S: 57-88-5
Form	Liquid
Storage Buffer	PBS, pH 7.4
Concentration	Lipid Concentration 10mM
Stability	2 months
Storage	2 °C to 8 °C (do not freeze)