



DOPC/CHOL Liposomes (100nm), DiA labeled;

Catalog: DOPC-014L

PRODUCT INFORMATION

Name	DOPC/CHOL Liposomes (100nm), DiA labeled;
Cat No.	DOPC-014L
Product Overview	The liposomes are prepared by the extrusion method to 100nm with very tight particle size distribution. The pr oduct is sterile filtered and filled in sterile vials guaranteed to be stable for 6 months when stored in refrigerato r. Lipid composition: DOPC/CHOL (54:45 mol/mol) Mean particle size: 100 nm (90-120 nm)
Lipid Composition	DOPC; CHOL
Application	Liposome production; Synthetic lipid
Storage Buffer	Hydration buffer: 10% sucrose, 20mM HEPES, pH 7.3 \pm 0.2
Concentration	Lipid concentration: 50 mM (50-55 mM) DiA: 0.5mM (0.39mg/mL)
Stability	6 Month for unopened vials.
Storage	Store at 2-8 centigrade.
Synonyms	DOPC; 1,2-dioleoyl-sn-glycero-3-phosphocholine; CHOL; cholesterol